

UG – Research Experiments

Module 1 - PCR

- Abhijeet Vikram

Overview

Objective - To Amplify a Gene with a SNP associated to a Human disease using PCR and its Extraction

Procedure -

- 1) Design best primers.
- 2) Find the Optimal annealing temp. for primers.
- 3) Perform a PCR on PCR to amplify the product.
- 4) Do Gel Extraction of the product and get the product concentration by nanodrop.
- 5) Send for Sequencing.

Telomerase Reverse Transcriptase – TERT Gene

- Telomerase reverse transcriptase (TERT) is the catalytic subunit of the **telomerase enzyme**, which plays a crucial role in maintaining the length of **telomeres** by addition of the telomere repeat TTAGGG.
- **Telomerase expression plays a role in cellular senescence**, as it is normally repressed in postnatal somatic cells resulting in progressive shortening of telomeres.
- **Deregulation of telomerase expression in somatic cells may be involved in oncogenesis.** Studies in mouse suggest that telomerase also participates in chromosomal repair, since de novo synthesis of telomere repeats may occur at double-stranded breaks. [provided by RefSeq, Jul 2008]

Primers Design

Criteria:

- Primer length ~ 15-30 bp
- G-C content ~ 40 – 60%
- Melting Temp. ~ 55 - 60°C
- Avoid Di-nucleotide repeats or single base runs could lead to slipping or hairpin loop structures.
- Primers shouldn't be complimentary to each other, could lead to formation of primer dimers and hairpin loop structures.

Type: Homozygous

Chromosome Position: 5:1287079

SNP Site: GAGACCCG[C/T]CTGGTGCA

SNP Accession No.: rs2853677

Forward Primer (5'→3'): AGACTTTGGGGCAATACAGGG

Reverse Primer (5'→3'): AACAGACCCCAATCCCCCA


Product Length: 289

PP	Template Strand	Length	Tm	GC%
Forward	Plus	21	59.99	52.38
Reverse	Minus	19	60.47	57.89

Method


- Prepared 20 μ L working stock of the primers using NFW.
- Prepared Master mix
 - Template DNA - 1.5 μ L (HEK Cells)
 - Forward Primer - 1 μ L
 - Reverse Primer - 1 μ L
 - DNTPs - 2 μ L
 - Taq polymerase - 1 μ L
 - 10 \times Buffer - 5 μ L
 - NFW – 38.5 μ L
- Performed PCR with 4 samples with a gradient of Annealing temp.

PCR – To find the optimal annealing temp.



Gel prep

- Prepared 1% agarose gel + 2 μ L EtBr (50ml)
- Took 8 μ L of 10 μ L from each sample and added 2 μ L of Loading Dye.
- Prepared Ladder (100bp) - 4 μ L Ladder + 4 μ L Loading Dye.
- Loaded the gel and ran at 120 V.






Fig 1. PCR results of all 4 samples.





Fig 2. Performed PCR on PCR on Sample 4

1. Performed PCR on Sample 4 again from template DNA (55°C) (35μL - Preserved)
2. Performed PCR on PCR on 2μL of sample 2 at 56.2°C (35μL - Preserved)



Gel Extraction

- Used 70 μ L (35 + 35 μ L) and ran it on 1% agarose gel at 80V and cut the product out.
- Used QIAquick Gel extraction kit to extract the product



Procedure:

- 1) Weighed the product
- 2) Add buffer QG and incubate at 50°C
- 3) Add isopropanol (100 μ L per 100mg)
- 4) Put it in QIAquick column and centrifuge for 1min
- 5) Discard flow through and wash with 0.75 ml of Buffer PE to QIAquick column and centrifuge for 1 min
- 6) Discard flow through and centrifuge for 1 min.
- 7) Elute DNA add 50 μ L of Buffer EB to the center of the QIAquick membrane and centrifuge the column for 1 min in a clean 1.5 ml microcentrifuge tube.

Nanodrop

- To measure the concentration of the result.
 - Select your desired application (DNA, RNA, or protein).
 - Gently clean and buff the pedestal with a lint-free wipe.
 - Pipette 1-2 μL of your blanking solution (Water in my case).
 - Click “Measure blank.”
 - Lift the arm and wipe clean with a lint-free wipe.
 - Then Put the 1-2 μL of sample and click sample.

- Results :

DNA concentration: 36.8 ng/ μL

A260 value: 0.737

A260/A280 ratio: 1.86

A260/A230 ratio: 0.26

DNA concentration: 16.7 ng/ μL

A260 value: 0.737

A260/A280 ratio: 1.81

A260/A230 ratio: 0.27

So, the final **DNA concentration of combined sample is 26.75 ng/ μL .**

Thank You!

Equipements used -

- NCBI Primer-BLAST – Primer designing
- HiMedia (100bp) – ladder
- QIAquick gel extraction kit - extraction
- Nanodrop